

# L&L SKIN


## MIO

- 2 -

Japan L&L Skin Tech Co. Limited  
No. 305, 2-24-8, Shimorenjaku, Mitaka, Tokyo, Japan, 181-0013


**\*Important Reminder: 1. Please charge to saturation before use.**  
**2. Pls recharge the device if it's not used for a long time(about 3 months).**

### Product Structure



### Product Features

- Design 680R & 700R massage arc line, fit face well, promote your cream and serum absorption, enjoy face massage & lifting.
- 3massage points, not only could fit face massage but also eye exerise, head massage. Even whole body massage. You could enjoy professional massage at home.
- IPX5 waterproof.
- Wireless charging, only need put on the customized charging base, easy-charge.
- Per 1minute intelligent automatic reminder function, easy to master the usage time.
- Portable, enjoy professional massage anytime, anywhere.



### Dual Modes For Selection

**Soothing mode:** Under blue light flash indication, cool your skin with a low vibration, soothe your skin.



**Heated mode:** Under red light flash indication, about 45°C warm your skin with a higher vibration, help your skin rejuvenation.

### Skin Care Tips

- Use together with your favorite essence, lotion, cream, eye cream, mask and other skin care products, enjoy skin rejuvenation.
- When you do face mask, you could use it promote essence absorption.

### How To Use

- Apply proper lotion to face.
- Long press on/off button, switch into soothing mode, get start skin care.
- Transfer swith off button into heated mode, deeply strengthen skin care effect.
- After use, use the paper towel clean or wash it by water, and dry it in a well-ventilated place.



### Charging Tips

- USB cable connect the charging base.
- Put device into charging bas directly, easy-charge.

### Cautions

- Before using the product, please read the manual carefully and keep it for future reference.
- The product has the function of heating, approximate temperature is in 45° C. Different people are in some difference when sensing the temperature, some are likely to feel hot.
- Portable, enjoy massage anywhere, anytime.
- When the battery is low, the red light which indicating the charging will automatically power off after flashing continuously for 30 seconds.
- Wireless charging, the LED light will be under blue indication flashes. 3-4 hours full of charge.
- Do not use the product while charging.
- If not used, please save it well after fully charged.
- Do not allow children to contact the product and accessories to avoid any danger.

### Q&A

The product suddenly in no vibration, does not work, what is the reason?

- Increase the dosage of the skin care products.
- Check whether the product electricity is in shortage or not.
- Check if the device starts correctly.
- Restart the product after turn off it.
- Please contact the brand distributor.

Heated Mode is not heating, what's the reason?

- Check whether the product electricity is in shortage or not.
- Check if the device starts correctly.
- Restart the product after turn off it.
- Please contact the brand distributor.

Can I use this product every day?

Yes, the vibration intensity and temperature adjustment of this device are intelligently designed by a professional team, it is suggested to use multi-times a day.

I am in allergic skin, can I use it?

This device specially uses the food grade stainless steel which in line with international standards, and food-grade silicone with high activity, it is proved to be safe, healthy and no any irritation to skin, the people in allergic skin please feel assured to use it every day.

How to charge?

USB cable connect charging bas, put the device into charging bas into charging condition.

How long can I use when fully charged?

Use ten minutes a day for two weeks.

Why can't I use the product while charging?

In order to make the product charge quickly to saturation, and it is not convenient to use it in handwhen connected to the USB charging cable. In order to give you a better product experience, it is hereby set.

Can I not clean the comb after use?

Unclean instrument, the surface is easy to breed bacteria, invade your skin, suggest be sure to clean and store, good use habits will prolong the service life of the product.

### Specification

Product name	MIO -2-
Voltage	5V
Rated Input Power	3W

### Warranty Card

▼ Fill In By Customer

Purchase Date \_\_\_\_\_

Name \_\_\_\_\_

Phone \_\_\_\_\_

Address \_\_\_\_\_

▼ Fill In By Seller

How To Buy ☐ Entity Shop ☐ Online Shop

Shop Name \_\_\_\_\_

Address \_\_\_\_\_

Contact Person \_\_\_\_\_ Phone \_\_\_\_\_

- The warranty period is one year, effectively from the date of purchase.
- The warranty scope only refers to the problems of L&L SKIN branded products caused due to the production, technique and quality.
- The cases of the following circumstances, even during the warranty period, should also be repaired by user cost
  - The damage caused by usage, maintaining and storage not according to user manual.
  - The damage caused due to incorrect operation, such as power is not normal or broken ( eg: drop)
  - The damage caused by dismantling of maintenance people who are not appointed by L&L SKIN brand.
  - Those who cannot present this warranty card and a valid purchase invoice;
  - Those Warranty recorded Information not conform to the actual commodity remark, or altered identification happens.
  - Damage due to force majeure (such as lightning, earthquakes, fire, flood, etc.)
- When contacting brand distributor or manufacturers for warranty repair, please present this warranty card and purchase invoice.
- During warranty period, the manufacturer will replace the faulty part free, but the freight out and home shall be borne by the customer.
- The warranty card shall be filled in block letters in clear written, no alteration allowed, please keep well, if lost, there will be no duplicate provided.
- The final interpretation of the warranty card belongs to Brand L&L SKIN.
- Due to the reason of Product upgrade, or package upgrade, the information contained in this warranty is subject to change without notice.